

INSECT BITE HEALER

STICK

Insect bite healer

NEW



2 Programmi
Settings

3 s
6 s

3 sec: for sensitive skin
6 sec: for normal skin



Gives relief



Without chemical
substances



2 Settings



5 year warranty

Medel **STICK** is a device that promotes healing of bites and stings caused by mosquitoes and insects, quickly soothing and reducing itching and swelling. The device features a ceramic plate that accelerates the healing.

- 2 application programs: 3 seconds on delicate skin and 6 seconds on normal skin
- Batteries included



Suitable for
pregnant women



Suitable for children
(3 years +)



Modern and
handy design

Activation button

ON-OFF button



TECHNICAL SPECIFICATIONS

Without chemical substances	√
Functions	2 time settings (3s/6s)
Max temperature	50°C ± 2°C (122°F ± 3,6°F)
Medical device	√
Power supply	2 x 1.5V AAA
Weight	Ca. 47g (with batteries)
Size	(LxWxH) 9.8 x 3.5 x 1.9 cm



6 pieces display box

95265

